

SixPACS Teleradiology Management System

User Guide Referring Physician Supplement

> Revision 2 04 Febuary 2014

	Table of Contents	
1. SixPACS Roles.	5	3
2. Life Cycle of ar	n Exam	3
3. Login		4
4. Search Exams.		4
5. Patient Detail		4
6. Viewing Image	es	4
7. Viewing Report	rts	5

SixPACS Roles

Access to the SixPACS system is organized into 8 different roles: Facility, Dispatcher (Staff), Tech, Radiologist, Referring Physician, Transcriptionist, Admin and Billing. Each role is capable of performing certain tasks, which follow the life cycle of an exam. Below are the definitions of each role.

Facility Role

The Facility user can view exams of their patients that have been read by a Radiologist, Transcribed or verified. This access is read only and they cannot change a record. The Facility user can also create new exams and track existing. They also receive the final report of the exam.

Dispatcher Role

The Dispatcher's role consists of several different tasks. They are responsible for fielding calls to schedule exams. If the patient doesn't have a detailed record in TMS, the dispatcher will create a patient record and then add the current exam. If the patient has a previous record and exam history, the Dispatcher will add the new exam to that patient's record.

Once the record is created for the patient and exam, the dispatcher is responsible for assigning the exam to a Tech. If the Tech is unable to perform the exam for any reason, the Dispatcher will need to re-assign the exam to another Tech. If an exam is canceled, the Dispatcher will mark the exam record as 'Canceled'. Dispatchers have the ability to move exam status backward using the 'Rollback' function, or forward by using the 'Exam Completed' function. Lastly, the Dispatcher is responsible for tracking the status of the exam and making sure it is completed in a timely fashion.

Tech Role

The Tech is responsible for uploading the images for new exams. The Tech can also verify the patient and exam info that was entered by dispatch and make any necessary changes. The Tech also has the ability to do add-on exams. They can search patient records, add a new patient or add an exam to an existing patient and have it assigned directly to them. This allows the Tech to operate independently after hours.

Radiologist Role

The Radiologist can use the TMS system to view the images attached to patient exams. They can view the images from a single exam, or they can also load images from past exams in order to do comparisons. After the images have been viewed, the Radiologist can mark them as 'Read' so that they are ready for the Transcriptionist. After the Transcriptionist has entered a transcription, the Radiologist can review it, revise it, and then verify it.

Referring Physician Role

The Referring Physician is allowed as an ordering physician to view completed exam images and reports, regardless of which facility the patient belongs. A user with the Physician role can view any patient that is "tagged" for them to access.

Transcriptionist Role

The Transcriptionist uses TMS to enter a typed transcription of the Radiologist's result and impression. They have the ability to revise previous transcriptions up until the point they have been verified by the Radiologist. The Transcriptionist can also view or print the transcription.

Admin Role

The Admin can be considered a manager who has access to all areas of the workflow of an exam. Their role is to make sure that all stages of an exam life cycle complete successfully. The Admin also has access to an Admin page where the supplemental tables such as Users, Facilities, Symptoms, etc. can be edited.

Billing Role

The Billing user has additional screens that are used to add billing information to the patient and exam records, to facilitate the billing process. This information includes selecting the correct procedure and symptom codes.

Life Cycle of an Exam in SixPACS

The life cycle of an Exam in SixPACS follows a path through as many as 13 different statuses. Depending on the user's role, they will be able to see exams at some or all of these stages.

Canceled - Canceled exams appear at the top of the worklist and will stay there for 24 hours. After 24 hours, they will be removed from the worklist.

To Be Scheduled – Certain exams, such as holter monitors, ultrasounds and pacer checks require additional scheduling, do to special equipment or patient preparation that must be done to perform the procedure. When these types of exams are ordered, they will be assigned this status, in order to allow staff to determine the date of service.

Call Facility - Once an exam has been scheduled in the "To Be Scheduled" status, it is changed to "Call Facility". In this status, the Staff user contacts the facility to verify that the date of service is acceptable. Also, any prep instructions for the patient are given at this time. Once the date of service is verified, the exam will change status to "Received", in order to be dispatched.

Received – Once a call has been received and the exam detail has been entered, the exam is moved to the 'Received' status and is ready to be assigned to a Tech.

Assigned – Exams are moved to 'Assigned' after they have been given to a tech. They will remain in the 'Assigned' status until the procedure is performed by the Tech. and all relevant information has to be received into the system.

Escalated – This status is used by techs when there is an issue transmission, or other problems arise which the tech needs assistance from management.

To Be Read – After the Tech has uploaded the image(s), the exam is automatically moved to the 'To Be Read' status and is ready to be reviewed by the Radiologist.

Order To Be Completed – Exam records that are read by outside radiology groups may have an external website that is utilized. This section shows exams that have had the images sent to the outside group, require an order be completed on the other side. Once the order is completed, the exam status is set to "To Be Transcribed".

To Be Transcribed – After a Radiologist has viewed the images associated with an exam and dictated the results and impression, it is marked and moved to the status of `To Be Transcribed'.

Q/A-

To Be Faxed – Once the Transcriptionist has typed in the Radiologist's impression and result, the exam status is changed to 'To Be Faxed'. If there is a fax location specified, the report will be automatically faxed to the facility.

To Be Called – Positive reports that are faxed are set to a status of "To Be Called". Once the user calls the facility to verify that they received the report, the user marks the exam as "Called". This sets the status of the exam to "To Be Billed".

To Be Billed – After the report has been faxed to the facility, the exam record is set to a status of "To Be Billed". These records will stay in the worklist for 24 hours, then they will be removed.

In SixPACS, the role of the Referring Physician provides access to view exam images and to review the transcribed report.

Login

Each Referring Physician will be able to login from any computer with a suitable MS Internet Explorer web browser. At the login screen, enter your username and password. The username will be your first initial and last name.

Search Exams

After logging in, the Referring Physician Menu will be displayed (Fig. 1). The Referring Physician Menu displays only 1 option: Search Exams

SixPACS Teleradiology Management System

Search Exams

Fig. 1 – Referring Physician Menu

Selecting the "Search Exams" menu option will display a search form to allow the Physician to search for exams that have been "tagged" for them to view (Fig. 2).

									SIXPACS Teleraulology Plan	lagement system
SEARCH										
										🖼 Worklist
SSN	PATIENT LAST	PATIENT FIRST	DOS FROM	DOS TO	EXAM ID	ORDER ID	PATIENT ID	Search		

Fig. 2 – Search Exams Page

Type any criteria into the search for to search for a particular exam, and click on the Search button. A list of exams will appear, as shown (Fig. 3). To view the details of a patient, click on the SSN. The Patient Detail will be shown (Fig. 4).

SEARCH											
											Search Worklist
SSN	PATIENT LAST	PATIENT FIRST	DOS FROM	DOS TO	EXAM ID	ORDER ID	PATIENT ID	Create			
Search Results) 23	Records Found							Search			
SSN (Click to view)	PATIENT LAST	PATIENT FIRST	DOB	FACILITY			TECH ASSIGNED	STATUS	DIAM ID	D05	REGION
303-34-1713	CAANFORD	1048	07131-27	0400108-040	CHAR-CENTER		788.41 (0.0.)484	100.0100	101420	10110-012	Simulation (
147 13 4775	TRANSPOLIS.	MADELINE	10/25/14	BLUE ADDRESS	CURNER BOUR	N LIVENDCROTHE	CHEV BALEMORE	16407080	104025	14/12/14	mentalisety
100-08-8779	040440-0.28	made.ind	0.01011-14	nes, for bits !!	a Admini		INCLU TROOM	viet of lash	201800	04/01/08	mantaliung
147-28-8772	DEPENDING 18	anddund	10/20/14	Manhood includes	OF HARRONSON		AND A SA MONDA	144.0100	100071	DOLLARD.	Maintakung 1
147-33-4712	100 anit-0.25	00006.018	100004	Manifelia and And	OF 1444108-041		Clinity 68, 57million	VAR INVEST	200400	10.000	Naniahorg.
107-28-8775	Demonstration of the	MADE IN	ADDRESS OF	Manhold Sectors	OF SAME DESIGNATION.		CLAT INLIGH	144.0100	275475	BUTCORF	manufactory
147-15-1712	Contaminut, 29	minfail, that	100.00114	Maritie within	OF HARRISON, AL		martial domains	van origin	2010.000	10.040	international design
VAR-DA APPR	Distantion, 18	water the	all the second	Maritie while	OF MARKING, ST.		AUTo 100164	LAN PROPERTY.	ansans!	RUNDER	International Contention of the
147-14-4111	Date and and	statisticity	10.00.14	witness with our	A # 2580.00 B			1488-001000	Long To	ALCONDA.	Reportation 10
and an interest	Tableson in the	ward do not	and the second	Lowing Print links	A REAL PROPERTY.			Land on Land	TRAMES.	-	manufactory 1
and the store	California, 14	martini dad	Andre Series and	Link hits as	A 4 10 10 10 10 10		Classics day, Sciences	when private	20100	ALC: NO.	distantion in the local distance of the loca
LAR DR ARRY	Concession in the	market and	and the second	Lota Analis	a local data		Date American	- and the second	and the second second	and the second	The second se
100.00.000	Concession in the	and a loss	and the second	And Alla and			and Chine	-			
100000000000000000000000000000000000000			10.00.10	Concernance and	and the second second second		State of State of State			10.04.03	
199.00.000				CONTRACT NO.	COLUMN ACCOUNTS		David Longian	-		ALC: NO.	action of the local division of the local di
181 11 122	1404 (180)	ALC: NOT THE OWNER OF THE OWNER OWNER OF THE OWNER OWNE OWNER OWNE OWNER OWNE OWNER OWN	100100-07	COUNTRY ME	solved stranged	IN STREET	BOUTT BOUTS	1000.000	Designed.	100.071.04	Activements

Fig. 3 – Search Results

								ps	Search
	554 187 05 9773	Patent Last		Patient First (Mobile 1, 1986)	008 3/20/19	18			
	Insun	ince		Number	Group		Auth Reg.		
	Primary MEDIO	LANS				51			
	Secondary Advis	A US HEALTHCARE		10100408305					
	Tertiary			-		- 10			
	Commente								
							View	Single	Exam
	-							Images	2
	Exam History		-	Click here for Comp	arisons			1	
	P View/Compare	Selected Exams	-	click here for comp	01150115		_	1	-
	05/07/2007	390690	VERIFE	ED				P (2)	6
1	Exam Type		Region		Auth #	Dispatcher	084		
/	Exerr(s) Perfor	med	Facility	na Mariana ana amin'ny fisiana	Physician	Tech Assig	ned .		T
/	CHEST IV		(458) 31	WISH HOME OF HARRIS	DR SAACKS	GELSINGE	L, CINDY	Mie	
/	Exam Requeste CHEST	нđ	105-PR	IVATE SOU WING	Requested Sy MATT	Radiologist RIGGID M.	D., CHARLES	vic.	IN INC
doxes	Symptoms CONCESTION		Call Tak	07 10:12:25 AM	Cenceled By				
exams	Comments		-	00000000000	1				
1									
-									
	03/15/2007	EXAM 10 380871	VERIFI	ED				P (2)	12
	Exam Type		Region		Auth #	Dispatcher		1	
	XRAY		Harrisbe	eng	-	Carr, Joyc			
	CHEST 1V	med.	(458) 38	WISH HOME OF HARRIS	SAACKS, DENNIS	MIXON, AN	IDREA		
	Exam Requests OREST	of.	Ndg./Room # 105 PRIVATE 5 WING Call Taken 3/12/2007 1/04/39 PM		Requested By	Radiologist SHITH D.O., HENRY			
	Bymptoma COMPARE 02/1	9/2007			Canceled By				
	Commanta								



View Exam Image

To view the images associated with the exam, click on the view icon (magnifying glass). A new window will open, the images will load into the viewer application. (refer to the Viewer Supplement for more information on the use of the viewer).

To view a comparison of the current exam with a prior exam, put a check in the checkbox next to the exam containing the images you would like to view. You can check more than one box if you would like to perform more than one comparison. Click on the "View/Compare Selected Exams" link at the top of the exam list. This will load the viewer in a new browser window. The image(s) will begin downloading immediately.

View Report

After the exam has been transcribed, the Physician has the ability to view the report. This is accomplished by clicking on the icon, as shown in Fig. 4. The report generated is in Adobe Acrobat, and requires the Adobe Reader to view the document. An example of a report is shown in Fig. 5.

				945 East Park Drive Suite 102 Harrisburg, PA 1711: Ph: (800) 420-9729 Fax: (717) 551-4467
PATIENT	Chearton			124 389 24 1753
FACILITY:	CRUE MILLER	NO CARE CEN	19.4	008 47/33/3987
ROOM:	45.0			
ORDERING PH	SICIAN	1 DE LARGE SCAR	*	
PROCEDURE:	KNEE LT 2V	KNEE RT 2V		
DOS:	03/15/200	7		
HISTORY:	UNSTABILI	Y, KNEE REPL	ACEMENTS	
RESULTS:				
In the left knee increased genu	the distal femur, pro there is likewise no l valgus deformity, po	cimal tobia, and cosening of the ssibly related to	posterior patella. In prosthetic componer instability.	ere is normal mild genu vaigu nts. However, there is slight
Diaba lassa Con		للاحد المحمد المح		
-	,			

Fig. 5 – View Transcription